

PRODUCT DATA SHEET
MEFIX®

Product name	Mefix®					
Product description	Self-adhesive fabric					
Intended use	Self-adhesive fabric for fixation of dressings, tubes and canulas etc.					
Product Codes	Art no	Description	Size (cm x m)	Size (inch)	Pack count con/ret/trp	Comments
	310250-05	Mefix 2,5 cm x 10 m	2,5 x 10	1 x 400	1/1/40	
	310276-05	Mefix 2,5 cm x 5 m	2,5 x 5	1 x 200	1/1/48	
	310299-00	Mefix 2,5 cm x 10 m, US	2,5 x 10	1 x 400	1/1/40	
	310500-05	Mefix 5 cm x 10 m	5 x 10	2 x 400	1/1/40	
	310540-01	Mefix 5 cm x 5 m, ph ES	5 x 5	2 x 200	1/1/24	
	310570-05	Mefix 5 cm x 2.5 m	5 x 2.5	2 x 100	1/1/30	
	310576-05	Mefix 5 cm x 5 m	5 x 5	2 x 200	1/1/30	
	310580-05	Mefix 5 cm x 2.5 m, ph FR	5 x 2.5	2 x 100	1/1/30	
	310584-05	Mefix 5 cm x 10 m, ph FR	5 x 10	2 x 400	1/1/40	
	310586-05	Mefix 5 cm x 5 m, ph FR	5 x 5	2 x 200	1/1/30	
	310599-00	Mefix 5 cm x 10 m, US	5 x 10	2 x 400	1/1/40	
	311000-05	Mefix 10 cm x 10 m	10 x 10	4 x 400	1/1/20	
	311040-01	Mefix 10 cm x 10 m, ph ES	10 x 10	4 x 400	1/1/30	
	311070-05	Mefix 10 cm x 2.5 m	10 x 2.5	4 x 100	1/1/30	
	311076-05	Mefix 10 cm x 5 m	10 x 5	4 x 200	1/1/30	
	311080-05	Mefix 10 cm x 2,5 m, ph FR	10 x 2,5	4 x 100	1/1/30	

	311084-05	Mefix 10 cm x 10 m, ph FR	10 x 10	4 x 400	1/1/20	
	311086-05	Mefix 10 cm x 5 m, ph FR	10 x 5	4 x 200	1/1/30	
	311090-05	Mefix 10 cm x 11 m, DE	10 x 11	4 x 440	1/1/20	
	311099-00	Mefix 10 cm x 10 m, US	10 x 10	4 x 400	1/1/20	
	311500-05	Mefix 15 cm x 10 m	15 x 10	6 x 400	1/1/22	
	311570-05	Mefix 2,5 m x 15 cm	15 x 2,5	6 x 100	1/1/20	
	311576-05	Mefix 15 cm x 5 m	15 x 5	6 x 200	1/1/20	
	311580-05	Mefix 2,5 m x 15 cm, ph FR	15 x 2,5	6 x 100	1/1/20	
	311586-05	Mefix 15 cm x 5 m, ph FR	15 x 5	6 x 200	1/1/20	
	311590-00	Mefix 15 cm x 11 m, DE	15 x 11	6 x 440	1/1/22	
	311599-00	Mefix 15 cm x 10 m, US	15 x 10	6 x 400	1/1/22	
	312000-05	Mefix 20 cm x 10 m	20 x 10	8 x 400	1/1/10	
	312076-05	Mefix 5 m x 20 cm	20 x 5	8 x 200	1/1/10	
	313000-05	Mefix 10 m x 30 cm	30 x 10	12 x 400	1/1/12	
	313076-05	Mefix 5 m x 30 cm	30 x 5	12 x 200	1/1/12	
Manufacturer	Mölnlycke Health Care AB Box 13080 SE-40252 Göteborg Sweden					
Country of Origin	Finland					
Product classification according to MDD	Class I					
CE-marked	Yes					
Sterilization method	N/A, Mefix [®] is a non sterile product					
Shelf life	5 years					
Product materials	Non woven	Spun laced non woven material that consist of polyester fibres				
	Adhesive	Water based poly acrylic adhesive				
	Release	Siliconised paper				
	Core (only for 310540-01 and 311040-01)	Polystyrene				

Other info of special interest	LATEX: Mefix [®] Products and packaging are latex free.			
Physical characteristics	Mefix [®] is released according to following standards:			
	Property	Value	Reference method	
	Water vapour permeability	>6500 g/m ² /24h	ASTM D 6701-01	
Packaging	The purpose of the dispenser box and transport box is to contain and protect material from dust, moisture and damage during transport, handling and storage.			
Storage Conditions	N/A			
GS1 Codes	Art No	CON	RET	TRP
	310250-05	-	7332551586225	7332551586218
	310276-05	-	7332551829247	7332551829230
	310299-00	-	7332551586263	7332551586256
	310500-05	-	7332551586294	7332551586287
	310540-01	-	7323190037186	7323190037179
	310570-05	-	7332551829605	7332551829599
	310576-05	-	7332551829681	7332551829674
	310580-05	-	7332551829643	7332551829636
	310584-05	-	7332551829797	7332551829780
	310586-05	-	7332551829759	7332551829742
	310599-00	-	7332551586379	7332551586362
	311000-05	-	7332551586409	7332551586393
	311040-01	-	7323190034888	7323190034871
	311070-05	-	7332551829827	7332551829810
	311076-05	-	7332551829889	7332551829872
	311080-05	-	7332551829858	7332551829841
	311084-05	-	7332551829940	7332551829933
	311086-05	-	7332551829919	7332551829902
	311090-05	-	7332551829971	7332551829964
	311099-00	-	7332551586447	7332551586430
	311500-05	-	7332551586485	7332551586478
	311570-05	-	7332551830007	7332551829995
	311576-05	-	7332551830069	7332551830052
	311580-05	-	7332551830038	7332551830021
	311586-05	-	7332551830090	7332551830083
	311590-00	-	7332430678546	7332430678539
311599-00	-	7332551586522	7332551586515	
312000-05	-	7332551586553	7332551586546	

	312076-05	-	7332551830120	7332551830113
	313000-05	-	7332551830151	7332551830144
	313076-05	-	7332551868567	7332551868550